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CELL TITER 96 AQUEOUS CELL PROLIFERATION ASSAY (PROMEGA)

1) Description :

The Cell Titer 96 Aqueous Cell proliferation assay is a colorimetric method for determining the number of viable cells in proliferation or chemosensitivity assay.

The Cell Titer 96 Aqueous assay is composed of solutions of a novel tetrazolium compound :

(3 (4,5-dimethylthiazol-2-yl) -5- (3-carboxymethoxyphenyl) 2- (4-sulfophenyl) -2H-tetrazolium, inner salt , **MTS**) and an electron coupling reagent (phenazine methosulfate ; **PMS**).

MTS is bioreduced by cells into a formazan that is soluble in tissue culture medium. The absorbance of the formazan at 490 nm can be measured directly from 96 well assay plates without additional processing. The conversion of **MTS** into the aqueous soluble formazan is accomplished by dehydrogenase enzymes found in metabolically active cells. The quantity of formazan product as measured by the amount of 490 nm absorbance is directly proportional to the number of living cells in culture.

2) Products :

One 96 wells assay plate,
PBS,

20 mL MTS solution (caution classified as irritants),

1 mL PMS solution (caution classified as a carcinogen),

Methyl hydroxybenzoate, propyl hydroxybenzoate solution 100 mg/mL

DMEM + 10 % FCS (foetal calf serum) medium,

Trypan blue 0,5 %,

Trypsin 0,25 % + EDTA solution,

and ELISA plate reader at 490 nm.

3) Protocol :

1) Maintain stock cultures of epithelial monolayer MS cells in DMEM medium containing 10 % of Foetal Calf Serum (FCS) in a 25 cm² flask.

2) Add 50 µL per well of hydroxybenzoate diluted in DMEM + 10 % FCS medium as described :

Quantity of hydroxybenzoate per well (µg)	Well
0	A1 → A4
20	B1 → B4
40	C1 → C4
60	D1 → D4
80	E1 → E4
100	F1 → F4
150	G1 → G4

3) Wash the MS cells with PBS. Add 1,5 mL Trypsin 0,25 % + EDTA solution. Incubated 1 min at 37°C to resuspend cells. Stop the enzyme action by adding 5 mL DMEM + 10 % FCS medium.

DOCUMENT N°6 : (suite)

- 4) Determine cells number and trypan blue viability and resuspend the cells to a final concentration of $1 \cdot 10^5 \cdot \text{mL}^{-1}$ in DMEM + 10 % FCS medium.
- 5) Dispense 50 μL of the cell suspension (5000 cells) into all wells of the plate prepared in step 2. The total volume in each well should now be 100 μL .
- 6) Incubate the plate for 48-72 hours at 37°C in humidified 5 % CO_2 atmosphere.
- 7) Add 20 μL per well of freshly prepared combined MTS/PMS solution (add 100 μL PMS solution to 2,0 mL MTS solution : gently swirl the tube to ensure complete mixing of combined MTS/PMS solution).
- 8) Incubate the plate for 1-4 hours at 37°C in a humidified 5 % CO_2 atmosphere.
- 9) Record the absorbance at 490 nm using an ELISA plate reader.

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