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Trypan Blue

Trypan blue is the most common stain used to distinguish viable cells from nonviable cells; only nonviable cells absorb the dye and appear blue and may also appear asymmetrical. Conversely, live, healthy cells appear round and refractile without absorbing the blue-colored dye. The use of this stain, however, is time- sensitive. Viable cells absorb Trypan blue over time, and can affect counting and viability results. Make dilutions just prior to counting to prevent viable cells from absorbing the stain, and thus appear non-viable.

In addition, trypan blue has a greater affinity for serum proteins than for cellular proteins. For cells cultured in high serum conditions, the background in the hemacytometer may be too dark. In this case, cells must be isolated from the serum. Simply centrifuge the cells and resuspend the cell pellet in a balanced salt solution or serum-free media prior to counting.

Procedure

1. Prepare a uniform cell suspension of the culture. Determine the number of cells per millileter and to be counted.

total number of cells using the following calculations:

2. Dilute a small sample of the cell suspension in 0.4% (w/v) Trypan blue. Other dilutions may be performed based on cell density of the suspension.

$$1:5 \text{ cells/ml} = \# \text{ of cells counted} \times 10^4 \times \text{dilution factor}$$

$$\# \text{ squares counted}$$

$$\text{total cells} = \text{cells/ml} \times \text{vol. of original cell suspension}$$

3. Center a cover glass over the hemacytometer. The percentage of viable cells can also be chambers. Fill one chamber with the cell dilution calculated using the following formula:

using a Pasteur pipette. The solution will pass under the cover glass by capillary action. Do not overfill. If the solution spreads into the two lateral grooves adjoining the grid table, clean the hemacytometer and repeat the application. If there are any bubbles in the solution covering the grid table, clean the hemacytometer and repeat the application.

$$\% \text{ Viability} = \frac{\# \text{ Viable Cells Counted}}{\text{Total \# Cells Counted}} \times 100$$

4. Place the hemacytometer on the stage of an inverted microscope and adjust focus using 100X magnification.

5. Use a hand-held counter to record cell counts in each of the four corner and central squares. Five squares (four corner and one center) are counted for a total of five squares. See Figure 1.

Important: Count squares touching the middle line of the triple line on the top and left of the squares. Do not count cells touching the middle line of the triple lines on the bottom or right side of the square.