Devoir n°3

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Document n°3 :

Detection of Enterotoxigenic Vibrio cholerae in Foods by the Polymerase Chain Reaction

U.S. Food & Drug Administration Center for Food Safety & Applied Nutrition 2001-OCT-24

Authors :Waiter H. Koch, William L. Payne, and Thomas A. Cebula

Recent epidemics of cholera in various parts of the world have emphasized the urgent need for rapid and reliable detection methods for *Vibrio cholerae*, especially in food and water. Classical microbiological methods are sensitive and specific; however, they require several days to complete and may result in considerable loss of perishable foods. Since cholera toxin production (encoded by the *ctxAB* genes) is the major factor in the pathogenesis of cholera, a polymerase chain reaction (PCR) method that selectively amplifies a DNA fragment within the *ctxAB* operon of *V. cholerae* has been developed and applied to various foods.

Many of these PCR methods have relied on extraction of DNA from contaminated foods, an additional step which adds several hours to the procedure and often requires modification of each diverse food matrix tested. An advantage of PCR is that the amplification reaction often proceeds well with crude lysates of cells, in some cases requiring only brief boiling of a bacterial suspension.

V. cholerae of the Inaba and Ogawa serotypes which lack the cholera toxin genes have been isolated; however, such strains are generally nonpathogenic. Since the presence of the cholera toxin operon is a prerequisite for pathogenicity, various PCR methods for the detection of pathogenic *V. cholerae* have all used the *ctxAB* genes as a target for amplification; these and the PCR method described here, will not detect nontoxic *V. cholerae*. As a practical matter, this PCR detection method allows one to define food samples as negative for the presence of toxic *V. cholerae* much more quickly than by following the complete microbiological identification scheme. However, it is recommended that alkaline peptone water (APW) enrichment broths used for PCR analysis

also be plated onto selective thiosulfate-citrate-bile salts-sucrose (TCBS) agar for isolation and direct confirmation of the presence of *V. cholerae* in samples that give positive PCR results.

A. Equipment and materials

1. For APW enrichment of *V. cholerae*.
2. Programmable automatic thermocycler
3. Horizontal gel electrophoresis apparatus
4. Electrophoresis constant-voltage powersupply
5. Heating plate
6. Microcentrifuge tubes, 1,5 and 0,6 mL
7. Variable digital micropipettors (e.g., 0,5-20 μ L, 20-200 μ L)
8. Aérosol-résistant pipet tips
9. Microcentrifuge
10. UV transilluminator
11. Polaroid camera
12. Polaroid film

B. Media and reagents

1. Alkaline peptone water (APW)
2. Cholera toxin PCR primers, 10pmol/ μ L stock solutions :
(5'-TGAAATAAAGCAGTCAGGTG-3'; 5'-GGTATTCTGCACACAAATCAG-3')
3. Taq DNA polymerase (native available from various vendors) 5 U/ μ L
4. 2'-deoxynucleoside-5'-triphosphates (dATP, dCTP, dGTP, dTTP); stock solution 1,25mM of each dNTP
5. 10X PCR buffer (100mM Tris-HCl, pH 8.3, 500mM KCl, 15mM $MgCl_2$)
6. Sterile deionized water
7. 10X TBE (0,9M Tris-borate, 0,02M EDTA, pH 8.3)
8. Agarose (nucleic acid electrophoresis grade)
9. Ethidium bromide solution, 10mg/ml
10. 6X sample loading buffer
11. DNA molecular weight markers (e.g., 123bp ladder, Bethesda Research Laboratories, Gaithersburg, MD)

C. Procedure for amplification of cholera toxin gene sequences from *V. cholerae* using APW enrichment broth

Food sample preparation and APW enrichment

APW enrichment lysate preparation.

Prepare APW washes or blends. Sample and freeze immediately (about one of four common buffers from Molecular Cloning: A Laboratory Manual 1mL) After enrichment (6-24 h), prepare crude APW lysates for PCR by Sambrook et al. (27) and load into sample wells of 1.5-1.8% agarose boiling 1mL samples in 1.5mL microcentrifuge tubes for approximately gel submerged in 1X TBE containing 1µg/ml ethidium bromide. After 5min. Lysates may be used for PCR immediately or stored in a -20°C appropriate migration with a constant voltage of 5-10 V/cm, illuminate the freezer until use.

NOTE: Due to the enormous amplification possible with the PCR, minute levels of contamination can result in false positives. It is recommended that sample preparation, PCR reaction set-up, and PCR product analysis be physically separated from one another to minimize contamination.

Minimally, use of aseptic technique in handling all PCR reagents and solutions is absolutely necessary. Use aerosol-resistant pipet tips for preparing samples and reagents for PCR reactions, and, if possible, a separate set of pipettors for analysis of PCR reaction products.

PCR reaction preparation.

To minimize cross contamination of PCR reagents, it is recommended that master mix solutions be prepared, aliquoted, and stored frozen. Master mixes contain all necessary reagents except Taq polymerase and that food. Minimally, this entails spiking 1 ml of a 1:10 and 1:100 APW food lysates being amplified. The final reaction contains 10mM Tris-HCl, pH 8.3; 50mM KCl, 1.5mM MgCl₂, 200mM each dATP, dCTP, dGTP, and dTTP; 2 equivalent amount of positive control lysate). A direct comparison of these to 5% (v/v) APW lysate; 0.5µM of each primer and 2.5U Taq polymerase spiked samples with the APW (no food) lysate containing identical per 100µL; reaction volumes of 25-100µL may be used. Add Taq polymerase to the master mix and add APW lysate upon distribution to 0.6mL microcentrifuge tube reaction vessels.

Temperature cycling.

While there is some variability in the heating and cooling dynamics of thermocyclers from different manufacturers, use of the following temperature cycling regimen should yield efficient amplification of the cfx gene fragment: Denaturation for 1min at 94°C, primer annealing for 1min at 55°C, and primer extension 72°C for 1min, repeated for no more than 35 cycles. Increasing the cycle number beyond 35 cycles often leads

Agarose gel analysis of PCR products.

Mix 10-20 pi portions of PCR reactions with 6X gel loading buffer (choose one of four common buffers from Molecular Cloning: A Laboratory Manual) and load into sample wells of 1.5-1.8% agarose gel submerged in 1X TBE containing 1µg/ml ethidium bromide. After appropriate migration with a constant voltage of 5-10 V/cm, illuminate the agarose gel with a UV transilluminator and visualize bands relative to molecular weight marker migration. The primers listed above give rise to a 777 bp fragment (19). Take Polaroid photographs of gels to document results. Further details regarding gel electrophoresis analyses may be found in the above-mentioned Molecular Cloning Laboratory Manual (27).

Proper controls.

The need for a number of control reactions to ensure accurate interpretation of PCR results cannot be overemphasized. Minimally, for PCR analysis of food types previously optimized for this method (e.g., vegetable washes, oyster, crab and shrimp blends), include a master mix contamination control containing no lysate and a toxigenic *V. cholerae* APW positive control in every analysis. For every new food blend to be analyzed by this PCR method, determine the potential inhibitory effects of that food. Minimally, this entails spiking 1 ml of a 1:10 and 1:100 APW food lysates being amplified. The final reaction contains 10mM Tris-HCl, pH 8.3; 50mM KCl, 1.5mM MgCl₂, 200mM each dATP, dCTP, dGTP, and dTTP; 2 equivalent amount of positive control lysate). A direct comparison of these to 5% (v/v) APW lysate; 0.5µM of each primer and 2.5U Taq polymerase spiked samples with the APW (no food) lysate containing identical per 100µL; reaction volumes of 25-100µL may be used. Add Taq polymerase to the master mix and add APW lysate upon distribution to either of the two food concentrations and prevents the occurrence of false negatives. It is unlikely that food washes (e.g., fruits and vegetables) will inhibit the PCR reaction unless the fruits are bruised and washing releases excessive acidity to the APW wash.

For additional information on this PCR method, contact Walter H. Koch at FDA, CFSAN, Division of Molecular Biological Research and Evaluation, 200 C St., S.W., Washington, DC 20204. Téléphone: (202) 205-4172 or (202) 205-5060; FAX: 205-4183; E-Mail: WHK@VAX8.CFSAN.FDA.GOV.